


Anti-Hsf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q00613

Item Number: CM0015415



Introduction

This rabbit polyclonal antibody targets HSF1 (HSTF1), validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Useful for stress response and transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	HSF1 Rabbit pAb
Remarks/Alias	HSTF1; HSF1
Species	Human
Gene ID	3297
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-529 of human HSF1 (NP_005517.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q00613
Protein Weight	57kDa
Shipping	Ice bag