

Anti-Hspa6 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P17066

Item Number: CM0028779



Introduction

HSPA6 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA; reacts with human, mouse, rat. Ideal for studying heat shock protein 70 family chaperone functions and stress response pathways.

[Learn More](#)

Parameter	Value
Product Name	HSPA6 Rabbit pAb
Remarks/Alias	HSP70B', HSPA6
Species	Human
Gene ID (Human)	3310
Gene ID	3310
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 474-643 of human HSPA6 (NP_002146.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17066
Protein Weight	71kDa
Shipping	Ice bag