


Anti-Human IgG Rabbit Monoclonal Antibody For Wb, Elisa - P01857 / P01859 / P01860 / P01861

Item Number: CM0003238



Introduction

Rabbit monoclonal antibody targeting human immunoglobulin heavy constant gamma chain (all subclasses). Validated for WB and ELISA, ideal for humoral immunity studies and diagnostic IgG quantification.

[Learn More](#)

Parameter	Value
Product Name	Human IgG Rabbit mAb
Species	Human
GeneID	3500
Immunogen	Synthetic Peptide Human IgG
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01857,P01859,P01860,P01861
Protein Weight	52 kDa
Shipping	Ice bag