


Anti-Human IL-5 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P05113

Item Number: CM0019541



Introduction

Rabbit polyclonal antibody against human IL-5, validated for WB, IF/ICC, ELISA. Cross-reacts with mouse, rat, and wide range predicted. Suitable for cytokine research, eosinophil studies, and B-cell differentiation assays.

[Learn More](#)

Parameter	Value
Product Name	Human IL-5 Rabbit pAb
Remarks/Alias	EDF, TRF, IL-5, Human IL-5
Species	Human
GeneID (Human)	3567
GeneID	3567
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-134 of human IL-5 (NP_000870.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P05113
Protein Weight	15kDa
Shipping	Ice bag