

Anti-Human Thioredoxin Reductase 1/Trxr1 Rabbit Monoclonal Antibody For Wb, Elisa - Q16881

Item Number: CM0005159



Introduction

Recombinant rabbit monoclonal antibody against human TXNRD1 (Thioredoxin reductase 1). Validated for Western blot and ELISA, cross-reacts with mouse and rat. Ideal for redox biology, cancer, and estrogen signaling research.

[Learn More](#)

Parameter	Value
Product Name	Human Thioredoxin Reductase 1/TRXR1 Rabbit mAb
Remarks/Alias	TR; TR1; TXNR; TRXR1; GRIM-12; Thioredoxin reductase 1 (TXNRD1)
Species	Human
Gene ID (Human)	7296
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-649 of human Thioredoxin reductase 1 (TXNRD1) (Q16881).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16881
Protein Weight	71kDa
Shipping	Ice bag