

Anti-Icam5 (Telencephalin) Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q60625

Item Number: CM0009046

Introduction

Rabbit monoclonal antibody targeting mouse ICAM5 (Telencephalin) validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Suitable for neuroscience research on neuron adhesion and dendritogenesis.

[Learn More](#)

Parameter	Value
Product Name	ICAM5 Rabbit mAb
Remarks/Alias	TLN; CD50; Tlcn; Icam3; ICAM-3; ICAM-5
Species	Mouse
GeneID (Human)	7087
GeneID	15898
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-227 of mouse ICAM5 (NP_032345.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q60625
Protein Weight	97kDa
Shipping	Ice bag