


Anti-Icos/Cd278 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y6W8

Item Number: CM0024340



Introduction

High-quality rabbit polyclonal antibody against human ICOS/CD278 (Q9Y6W8), validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for T-cell costimulation and immunology research.

[Learn More](#)

Parameter	Value
Product Name	ICOS/CD278 Rabbit pAb
Remarks/Alias	AILIM; CD278; CVID1; ICOS/CD278
Species	Human
GeneID (human)	29851
GeneID	29851
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 21-140 of human ICOS/CD278 (NP_036224.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6W8
Protein Weight	23kDa
Shipping	Ice bag