

Anti-Idh3A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P50213

Item Number: CM0018002



Introduction

High-quality rabbit polyclonal antibody against human IDH3A, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Recombinant protein immunogen ensures specificity. Suitable for studying TCA cycle and related diseases.

[Learn More](#)

Parameter	Value
Product Name	IDH3A Rabbit pAb
Remarks/Alias	RP90; IDH3A
Species	Human
GeneID	3419
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-366 of human IDH3A (NP_005521.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50213
Protein Weight	40kDa
Shipping	Ice bag