


Anti-Ifnar1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P17181

Item Number: CM0023341



Introduction

High-quality rabbit polyclonal antibody targeting human IFNAR1, validated for WB, IHC-P, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for type I interferon signaling research.

[Learn More](#)

Parameter	Details
Product Name	IFNAR1 Rabbit pAb
Remarks/Alias	AVP; IFRC; IFNAR; IFNBR; IMD106; IFN-alpha-REC; IFNAR1
Species	Human
Gene ID	3454
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 470-557 of human IFNAR1 (NP_000620.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17181
Protein Weight	64kDa
Shipping	Ice bag