


Anti-Igf2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P01344

Item Number: CM0021502



Introduction

IGF2 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, and ELISA applications with cross-reactivity for human, mouse, and rat samples. Recombinant immunogen ensures reliable detection of insulin-like growth factor II in fetal growth, metabolism, and differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	IGF2 Rabbit pAb
Remarks/Alias	GRDF; SRS3; IGF-II; PP9974; C11orf43; IGF2
Species	Human
GeneID (Human)	3481
GeneID	3481
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 25-180 of human IGF2 (NP_000603.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01344
Protein Weight	20kDa
Shipping	Ice bag