


Anti-Igf2Bp3 Polyclonal Antibody For Wb, Elisa - 000425

Item Number: CM0028782



Introduction

IGF2BP3 Rabbit pAb for WB, ELISA. Validated in Human, Mouse, Rat. Targets insulin-like growth factor 2 mRNA-binding protein 3 (IMP3/KOC). Ideal for cancer research and RNA biology studies.

[Learn More](#)

Parameter	Value
Product Name	IGF2BP3 Rabbit pAb
Remarks/Alias	KOC, CT98, IMP3, KOC1, IMP-3, VICKZ3
Species	Human
GeneID	10643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human IGF2BP3 (NP_006538.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O00425
Protein Weight	64kDa
Shipping	Ice bag