

Anti-Ikk α Rabbit Monoclonal Antibody For Wb, Elisa - O15111

Item Number: CM0007306



Introduction

Rabbit monoclonal antibody against human IKK α (CHUK, IKK1, IKKA) for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Ideal for studying NF-kappa-B signaling in immune and inflammatory response pathways.

[Learn More](#)

Parameter	Value
Product Name	IKK α Rabbit mAb
Remarks/Alias	BPS2; IKK1; IKKA; IKBKA; TCF16; NFKB1A; IKK-alpha; IKK α
Species	Human
Gene ID	1147
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O15111
Protein Weight	85kDa
Shipping	Ice bag