

Anti-IL-17Rb Monoclonal Antibody For Fc, Elisa - Q9Nrm6

Item Number: CM0006397



Introduction

Rabbit monoclonal antibody targeting human IL-17RB for FC and ELISA. Ideal for studying pro-inflammatory cytokine signaling, hematopoietic cell differentiation, and endocrine tissue biology. - Q9NRM6

[Learn More](#)

Parameter	Value
Product Name	IL-17RB Rabbit mAb
Remarks/Alias	CRL4, EVI27, IL17BR, IL17RH1
Species	Human
Gene ID	55540
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-289 of human IL-17RB (NP_061195.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC, ELISA
Cross Reactivity	Human
SWISS	Q9NRM6
Protein Weight	56kDa
Shipping	Ice bag