

Anti-IL-31 Monoclonal Antibody For Wb, If/Icc, Elisa - Q6Ebc2

Item Number: CM0007125

Introduction

Rabbit monoclonal antibody against human IL-31, validated for WB, IF/ICC, ELISA.

Recombinant immunogen, 18 kDa. Ideal for cytokine signaling and skin immunity research.

[Learn More](#)

Parameter	Value
Product Name	IL-31 Rabbit mAb
Remarks/Alias	IL-31
Species	Human
Gene ID (Human)	386653
Gene ID	386653
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-164 of human IL-31 (NP_001014358.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q6EBC2
Protein Weight	18kDa
Shipping	Ice bag