


Anti-IL11Ra Rabbit Polyclonal Antibody For Wb And Elisa - Q14626

Item Number: CM0028180



Introduction

IL11RA Rabbit Polyclonal Antibody validated for Western blotting and ELISA. Detects human IL11RA with cross-reactivity to mouse and rat. Suitable for research into IL11 signaling, craniofacial development, and cancer biology.

[Learn More](#)

Parameter	Value
Product Name	IL11RA Rabbit pAb
Remarks/Alias	CRSDA; IL11RA
Species	Human
Gene ID (Human)	3590
Gene ID	3590
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-290 of human IL11RA (NP_001136256.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14626
Protein Weight	45kDa
Shipping	Ice bag