

Anti-IL1Rn Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P25085

Item Number: CM0019256

Introduction

Rabbit polyclonal antibody against mouse IL1RN, an anti-inflammatory interleukin-1 receptor antagonist. Suitable for WB, IF/ICC, ELISA with cross-reactivity to human and mouse IL1RN. Ideal for studying immune regulation and inflammation.

[Learn More](#)

Parameter	Value
Product Name	IL1RN Rabbit pAb
Remarks/Alias	IL-1ra; F630041P17Rik; IL1RN
Species	Mouse
GeneID	16181
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 27-178 of mouse IL1RN (NP_001034790.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P25085
Protein Weight	20kDa
Shipping	Ice bag