

Anti-IL1B Polyclonal Antibody For Wb, Elisa - P01584

Item Number: CM0023312



Introduction

High-quality rabbit polyclonal antibody against human IL-1 beta, a key pro-inflammatory cytokine. Validated for WB and ELISA with cross-reactivity to mouse. Ideal for inflammation research.

[Learn More](#)

Parameter	Value
Product Name	IL1 β Rabbit pAb
Remarks/Alias	IL-1; IL1F2; IL1beta; IL1-BETA; IL1 β
Species	Human
GeneID	3553
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-170 of human IL1 β (NP_000567.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01584
Protein Weight	31kDa
Shipping	Ice bag