

Anti-IL21 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Hbe4

Item Number: CM0005678



Introduction

IL21 Rabbit Monoclonal Antibody validated for WB and ELISA applications. Reacts with human IL21, a key cytokine in Tfh cell generation and antibody responses. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	IL21 Rabbit mAb
Remarks/Alias	Za11; IL-21; CVID11; IL21
Species	Human
Gene ID (Human)	59067
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-162 of human IL21 (NP_068575.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HBE4
Protein Weight	17KDa/18KDa
Shipping	Ice bag