

Anti-IL7 Rabbit Polyclonal Antibody For Wb, Elisa - P13232

Item Number: CM0022153



Introduction

Rabbit polyclonal antibody targeting human Interleukin-7 (IL-7), validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studies on lymphocyte development, T-cell and B-cell homeostasis, and cytokine signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	IL7 Rabbit pAb
Remarks/Alias	IL-7; IL7
Species	Human
GeneID (Human)	3574
GeneID	3574
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 26-177 of human IL7 (NP_000871.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P13232
Protein Weight	20kDa
Shipping	Ice bag