

Anti-ILf3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q12906

Item Number: CM0021280



Introduction

Rabbit polyclonal antibody targeting ILF3 (NF90), validated for Western blot, IHC-P, IF/ICC, IP, and ELISA in human, mouse, and rat species. Ideal for studying RNA-binding proteins in circRNA biogenesis, transcriptional regulation, and antiviral defense.

[Learn More](#)

Parameter	Value
Product Name	ILF3 Rabbit pAb
Remarks/Alias	CBTF; DRBF; MMP4; MPP4; NF90; NFAR; NF110; NF90a; NF90b; NF90c; NFAR2; TCP80; DRBP76; NF110b; NFAR-1; NFAR-2; NFAR90; TCP110; MPP4110; NF90ctv; NFAR110; MPHOSPH4; NF-AT-90; ILF3
Species	Human
GeneID (human)	3609
GeneID	3609
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human ILF3 (NP_036350.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12906
Protein Weight	95kDa
Shipping	Ice bag