

Anti-Inhibin Alpha (Inha) Rabbit Pab For Wb, Ihc-P, Elisa - P05111

Item Number: CM0024372

Introduction

Rabbit polyclonal antibody targeting human Inhibin alpha (INHA), validated for WB, IHC-P, ELISA, with cross-reactivity to mouse and rat. Essential for reproductive endocrinology and HPG axis research.

[Learn More](#)

Parameter	Value
Product Name	Inhibin alpha (INHA) Rabbit pAb
Remarks/Alias	INHA; Inhibin alpha (INHA)
Species	Human
GeneID (Human)	3623
GeneID	3623
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human Inhibin alpha (INHA) (NP_002182.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05111
Protein Weight	40kDa
Shipping	Ice bag