


Anti-Inhibin Beta A (Inhba) Monoclonal Antibody For Wb, If/Icc, Elisa - P08476

Item Number: CM0012417



Introduction

Rabbit monoclonal antibody targeting human Inhibin beta A (INHBA) for WB, IF/ICC, ELISA applications. Cross-reacts with mouse and rat. Supports research in reproductive endocrinology, stem cell differentiation, and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	Inhibin beta A (INHBA) Rabbit mAb
Remarks/Alias	EDF; FRP; Inhibin beta A (INHBA)
Species	Human
Gene ID	3624 (Human)
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human INHBA (P08476)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08476
Protein Weight	47 kDa
Shipping	Ice bag