

Anti-Inos Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P29477

Item Number: CM0012551



Introduction

High-quality Anti-iNOS Rabbit Monoclonal Antibody validated for WB, IF/ICC, and ELISA. Reacts with mouse, human, and rat iNOS (Swiss-Prot P29477). Ideal for studying nitric oxide signaling and inflammation.

[Learn More](#)

Parameter	Value
Product Name	iNOS Rabbit mAb
Remarks/Alias	iNOS; Nos-2; Nos2a; i-NOS; NOS-II; MAC-NOS
Species	Mouse
Gene ID (Human)	4843
Gene ID	18126
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of mouse iNOS (NP_035057.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29477
Protein Weight	131 kDa
Shipping	Ice bag