

Anti-Integrin Alpha 5 (Itga5/Cd49E) Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P08648

Item Number: CM0028915



Introduction

Rabbit polyclonal antibody against human Integrin alpha 5 (ITGA5/CD49e). Validated for WB, IF-P, ELISA. Reacts with human, mouse, rat. Ideal for studying cell adhesion, integrin signaling, and fibronectin receptor functions.

[Learn More](#)

Parameter	Value
Product Name	Integrin alpha 5 (ITGA5/CD49e) Rabbit pAb
Remarks/Alias	FNRA; CD49e; VLA-5; VLA5A; Integrin alpha 5 (ITGA5/CD49e)
Species	Human
GeneID (Human)	3678
GeneID	3678
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 750-1000 of human Integrin alpha 5 (ITGA5/CD49e) (NP_002196.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08648
Protein Weight	115kDa
Shipping	Ice bag