

## Anti-Integrin-B1/Cd29 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Fc, Elisa - P05556

Item Number: CM0006621



### Introduction

Rabbit monoclonal antibody against human Integrin beta-1 (CD29). Validated for WB, IHC-P, IF/ICC, FC, and ELISA; cross-reacts with mouse and rat. Ideal for cell adhesion and signaling research.

[Learn More](#)

| Parameter        | Value                                                                                                                                            |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Integrin-β1/CD29 Rabbit mAb                                                                                                                      |
| Remarks/Alias    | CD29; FNRB; MDF2; VLAB; GPIIA; MSK12; VLA-BETA; Integrin-β1/CD29                                                                                 |
| Species          | Human                                                                                                                                            |
| GeneID (Human)   | 3688                                                                                                                                             |
| GeneID           | 3688                                                                                                                                             |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-728 of human Integrin-β1/CD29 (NP_596867.1). |
| Source           | Rabbit                                                                                                                                           |
| Category         | Monoclonal Antibodies                                                                                                                            |
| Application      | WB, IHC-P, IF/ICC, FC, ELISA                                                                                                                     |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                                |
| SWISS            | P05556                                                                                                                                           |
| Protein Weight   | 88 kDa                                                                                                                                           |
| Shipping         | Ice bag                                                                                                                                          |