


Anti-Integrin B5/Itgb5 Rabbit Monoclonal Antibody For Wb And Elisa - P18084

Item Number: CM0006799



Introduction

High-quality rabbit monoclonal antibody against human Integrin β 5 (ITGB5), validated for Western blot and ELISA. Recombinant immunogen (aa 24-490). Ideal for cell adhesion and viral entry research - P18084.

[Learn More](#)

Parameter	Value
Product Name	Integrin β 5/ITGB5 Rabbit mAb
Remarks/Alias	ITGB5; integrin beta-5; Integrin β 5/ITGB5
Species	Human
GeneID	3693
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-490 of human Integrin β 5/ITGB5 (NP_002204.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P18084
Protein Weight	88kDa
Shipping	Ice bag