

Anti-Interferon Alpha 1 (Ifn-A1) Monoclonal Antibody For Wb, Elisa - P01562

Item Number: CM0003407

Introduction

Rabbit monoclonal antibody against human Interferon alpha 1 (IFN- α 1), validated for Western blot and ELISA. Ideal for studying interferon signaling and immune response research.

[Learn More](#)

Parameter	Value
Product Name	Interferon alpha 1 (IFN- α 1) Rabbit mAb
Remarks/Alias	IFL; IFN; IFNA@; IFNA13; IeIF D; IFN-ALPHA; IFN-alphaD; Interferon alpha 1 (IFN- α 1)
Species	Human
GeneID	3439 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-189 of human Interferon alpha 1 (IFN- α 1) (NP_076918.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01562
Protein Weight	22kDa
Shipping	Ice bag