

Anti-Interferon Alpha 1 (Ifn-A1) Rabbit Polyclonal Antibody For Wb, Elisa - P01562

Item Number: CM0015455

Introduction

Rabbit polyclonal antibody against human Interferon alpha 1 (IFN- α 1), suitable for WB and ELISA. Recombinant immunogen corresponding to amino acids 24-189, validated for human samples. Ideal for cytokine research.

[Learn More](#)

Parameter	Value
Product Name	Interferon alpha 1 (IFN- α 1) Rabbit pAb
Remarks/Alias	IFL; IFN; IFNA@; IFNA13; IeIF D; IFN-ALPHA; IFN-alphaD; Interferon alpha 1 (IFN- α 1)
Species	Human
Gene ID (Human)	3439
Gene ID	3439
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-189 of human Interferon alpha 1 (IFN- α 1) (NP_076918.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01562
Protein Weight	22kDa
Shipping	Ice bag