


Anti-Interferon Alpha-2 (Ifn-A2) Monoclonal Antibody For Wb, Ihc-P, Elisa - P01563

Item Number: CM0003115



Introduction

Rabbit monoclonal anti-Interferon alpha-2 (IFN- α 2) antibody, validated for WB, IHC-P, ELISA in human, mouse, rat. Recombinant immunogen, 22kDa, SWISS P01563.

[Learn More](#)

Parameter	Value
Product Name	Interferon alpha 2 (IFN- α 2) Rabbit mAb
Remarks/Alias	IFNA; IFNA2B; IeIF A; IFN-alphaA; IFN-alpha-2; Interferon alpha 2 (IFN- α 2)
Species	Human
Gene ID (Human)	3440
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-188 of human Interferon alpha 2 (IFN- α 2) (NP_000596.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01563
Protein Weight	22kDa
Shipping	Ice bag