

Anti-Ipo4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Tex9

Item Number: CM0018004



Introduction

Rabbit polyclonal antibody against human Importin-4 (IPO4), suitable for WB, IF/ICC, and ELISA. Recognizes human, mouse, and rat IPO4 (~119 kDa). Ideal for nuclear transport and histone import research.

[Learn More](#)

Parameter	Value
Product Name	IPO4 Rabbit pAb
Remarks/Alias	Imp4; IPO4
Species	Human
Gene ID (Human)	79711
Gene ID	79711
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 728-1081 of human IPO4 (NP_078934.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TEX9
Protein Weight	119kDa
Shipping	Ice bag