

Anti-Irf3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P70671

Item Number: CM0019359



Introduction

Rabbit polyclonal antibody targeting mouse Irf3 (IRF-3, IIAE7), validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for innate immunity and interferon signaling studies.

[Learn More](#)

Parameter	Value
Product Name	Irf3 Rabbit pAb
Remarks/Alias	IRF-3, IIAE7, IRF 3
Species	Mouse
Gene ID (Human)	3661
Gene ID	54131
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 100-200 of mouse Irf3(NP_058545.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P70671
Protein Weight	46kDa
Shipping	Ice bag