

Anti-Irs2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y4H2

Item Number: CM0029071



Introduction

Rabbit polyclonal antibody targeting human IRS2, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Essential for insulin signaling, glucose metabolism, and cell proliferation research.

[Learn More](#)

Parameter	Value
Product Name	IRS2 Rabbit pAb
Remarks/Alias	IRS-2; IRS2
Species	Human
GeneID (human)	8660
GeneID	8660
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1238-1338 of human IRS2 (NP_003740.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y4H2
Protein Weight	137kDa
Shipping	Ice bag