

Anti-Itga3 Polyclonal Antibody For Wb, Elisa - P26006

Item Number: CM0019322

Introduction

Rabbit polyclonal antibody against human ITGA3, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for integrin alpha-3 function, cell adhesion, and cancer invasion research.

[Learn More](#)

Parameter	Value
Product Name	ITGA3 Rabbit pAb
Remarks/Alias	JEB7; VL3A; CD49C; FRP-2; GAPB3; ILNEB; MSK18; VCA-2; VLA3a; GAP-B3; ITGA3
Species	Human
GeneID (Human)	3675
GeneID	3675
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 85-260 of human ITGA3 (P26006).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P26006
Protein Weight	117kDa
Shipping	Ice bag