

Anti-Jam2 Rabbit Polyclonal Antibody For Wb And Elisa - P57087

Item Number: CM0026755

Introduction

Anti-JAM2 rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Cross-reacts with human, mouse, and rat JAM2. Ideal for studying leukocyte extravasation, hematopoietic stem cell homing, and tight junction biology.

[Learn More](#)

Parameter	Value
Product Name	JAM2 Rabbit pAb
Remarks/Alias	JAMB; CD322; IBGC8; JAM-B; VEJAM; PRO245; VE-JAM; C21orf43; JAM2
Species	Human
Gene ID (Human)	58494
Gene ID	58494
Immunogen	Recombinant fusion protein containing amino acids 29-238 of human JAM2 (NP_067042.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P57087
Protein Weight	33kDa
Shipping	Ice bag