

Anti-Jarid2 Rabbit Polyclonal Antibody For Wb And Elisa - Q92833

Item Number: CM0030707



Introduction

Rabbit polyclonal antibody against JARID2, validated for Western blot (WB) and ELISA. Cross-reacts with human, mouse, and rat. Targets a key PRC2 regulator involved in embryonic development and chromatin modification. Ideal for epigenetics research.

[Learn More](#)

Parameter	Value
Product Name	JARID2 Rabbit pAb
Remarks/Alias	JMJ; DIDDF; JARID2
Species	Human
Gene ID (Human)	3720
Gene ID	3720
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1027-1246 of human JARID2 (NP_004964.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92833
Protein Weight	139kDa
Shipping	Ice bag