


Anti-Jmjd1C Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15652

Item Number: CM0023332



Introduction

High-quality rabbit anti-JMJD1C polyclonal antibody validated for WB, IHC-P, ELISA. Recognizes human, mouse, rat JMJD1C (KDM3C/TRIP8), a lysine demethylase involved in DNA repair and transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	JMJD1C Rabbit pAb
Remarks/Alias	KDM3C; TRIP8; TRIP-8; JMJD1C
Species	Human
Gene ID (Human)	221037
Gene ID	221037
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 354-448 of human JMJD1C (NP_116165.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15652
Protein Weight	285kDa
Shipping	Ice bag