


Anti-Jnk3 Polyclonal Antibody For Wb And Elisa - P53779

Item Number: CM0021124



Introduction

Rabbit polyclonal antibody against human JNK3 (MAPK10), validated for Western blot and ELISA, with cross-reactivity to mouse. Ideal for studying stress-activated kinase signaling, neuronal apoptosis, and circadian rhythm regulation.

[Learn More](#)

Parameter	Value
Product Name	JNK3 Rabbit pAb
Remarks/Alias	JNK3; JNK3A; PRKM10; SAPK1b; p493F12; p54bSAPK
Species	Human
Gene ID (Human)	5602
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P53779
Protein Weight	44kDa
Shipping	Ice bag