


Anti-Jnk3 Rabbit Polyclonal Antibody For Wb, Elisa - P53779

Item Number: CM0022544



Introduction

Rabbit polyclonal antibody against human JNK3 (MAPK10), validated for Western blot and ELISA. Detects endogenous JNK3 in human, mouse, and rat samples. Useful for studying neuronal apoptosis and stress signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	JNK3 Rabbit pAb
Remarks/Alias	JNK3; JNK3A; PRKM10; SAPK1b; p493F12; p54bSAPK
Species	Human
Gene ID (Human)	5602
Gene ID	5602
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 300-400 of human JNK3 (NP_620448.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53779
Protein Weight	44kDa
Shipping	Ice bag