

Anti-Josd1 Polyclonal Antibody For Wb, Elisa - Q15040

Item Number: CM0013935

Introduction

Rabbit polyclonal antibody targeting human JOSD1 (Josephin-1) for WB and ELISA applications. Cross-reacts with mouse and rat. Useful for studying deubiquitination and membrane dynamics. Ideal for research on protein degradation and cell motility.

[Learn More](#)

Parameter	Value
Product Name	JOSD1 Rabbit pAb
Remarks/Alias	dj508l15.2; JOSD1
Species	Human
GeneID	9929
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15040
Protein Weight	23kDa
Shipping	Ice bag