


Anti-JunD Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P17535

Item Number: CM0018683



Introduction

Rabbit polyclonal antibody specific to JunD (AP-1 subunit) validated for WB, IHC-P, ELISA. Detects human, mouse, and rat JunD. Ideal for transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	JunD Rabbit pAb
Remarks/Alias	AP-1; JunD
Species	Human
Gene ID (Human)	3727
Gene ID	3727
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human JunD (NP_005345.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17535
Protein Weight	35kDa
Shipping	Ice bag