

Anti-Kaiso/Zbtb33 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q86T24

Item Number: CM0016497



Introduction

Rabbit polyclonal antibody against human Kaiso/ZBTB33, validated for WB, IHC-P, ELISA with cross-reactivity to mouse. Ideal for research on transcriptional regulation, Wnt signaling, and DNA methylation.

[Learn More](#)

Parameter	Value
Product Name	Kaiso/ZBTB33 Rabbit pAb
Remarks/Alias	ZNF348; ZNF-kaiso; Kaiso/ZBTB33
Species	Human
GeneID (Human)	10009
GeneID	10009
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 523-672 of human Kaiso/ZBTB33 (NP_006768.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86T24
Protein Weight	74kDa
Shipping	Ice bag