

Anti-Kat2A (Gcn5) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q92830

Item Number: CM0020348



Introduction

KAT2A rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat KAT2A (GCN5), a histone acetyltransferase involved in epigenetic regulation, transcription activation, and stem cell pluripotency. Ideal for chromatin research.

[Learn More](#)

Parameter	Value
Product Name	KAT2A Rabbit pAb
Remarks/Alias	GCN5; hGCN5; GCN5L2; PCAF-b; KAT2A
Species	Human
GeneID (Human)	2648
GeneID	2648
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92830
Protein Weight	94kDa
Shipping	Ice bag