


Anti-Kcnd3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uk17

Item Number: CM0024809



Introduction

Rabbit polyclonal antibody targeting human KCND3 (KV4.3), suitable for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat samples. Ideal for studying A-type potassium channel function in brain and heart.

[Learn More](#)

Parameter	Value
Product Name	KCND3 Rabbit pAb
Remarks/Alias	KV4.3; SCA19; SCA22; BRGDA9; KCND3L; KCND3S; KSHIVB; KCND3
Species	Human
GeneID (Human)	3752
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 521-655 of human KCND3 (NP_004971.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UK17
Protein Weight	73kDa
Shipping	Ice bag