

Anti-Kcnj6 Rabbit Polyclonal Antibody For Wb And Elisa - P48051

Item Number: CM0029276

Introduction

Rabbit polyclonal antibody targeting human KCNJ6/GIRK2 inward rectifier potassium channel. Validated for WB and ELISA, cross-reacts with human and mouse. Suitable for insulin secretion and ion channel studies.

[Learn More](#)

Parameter	Value
Product Name	KCNJ6 Rabbit pAb
Remarks/Alias	BIR1; GIRK2; KATP2; KCNJ7; KPLBS; GIRK-2; KATP-2; KIR3.2; hiGIRK2; KCNJ6
Species	Human
Gene ID (Human)	3763
Gene ID	3763
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human KCNJ6 (NP_002231.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48051
Protein Weight	48kDa
Shipping	Ice bag