

Anti-Kdm3A Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Y4C1

Item Number: CM0013304

Introduction

This rabbit polyclonal antibody detects human KDM3A (JMJD1A, TSGA), a 147kDa histone H3 Lys-9 demethylase. Suitable for WB, IHC-P, IF/ICC, and ELISA. Ideal for chromatin and transcription research.

[Learn More](#)

Parameter	Value
Product Name	KDM3A Rabbit pAb
Remarks/Alias	TSGA; JMJD1; JHDM2A; JHMD2A; JMJD1A; KDM3A
Species	Human
GeneID (Human)	55818
GeneID	55818
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 975-1283 of human KDM3A (NP_060903.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y4C1
Protein Weight	147kDa
Shipping	Ice bag