


Anti-Kdm5C Rabbit Polyclonal Antibody For Wb And Elisa - P41229

Item Number: CM0013097



Introduction

High-quality rabbit polyclonal antibody against human KDM5C, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying histone demethylation and transcriptional regulation. SWISS: P41229.

[Learn More](#)

Parameter	Value
Product Name	KDM5C Rabbit pAb
Remarks/Alias	MRXJ; SMCX; MRX13; MRXSJ; XE169; MRXSJ; JARID1C; DXS1272E; KDM5C
Species	Human
Gene ID (Human)	8242
Gene ID	8242
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 1400 to the C-terminus of human KDM5C (XP_011529128.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P41229
Protein Weight	176kDa
Shipping	Ice bag