


Anti-Ki67 Polyclonal Antibody For Wb, Ihc-P, Elisa - P46013

Item Number: CM0023343



Introduction

Rabbit polyclonal antibody against human Ki67 (proliferation marker protein). Validated for western blot, immunohistochemistry (paraffin), and ELISA. Cross-reacts with mouse. Designed to support cell proliferation assays and cancer biomarker discovery.

[Learn More](#)

Parameter	Value
Product Name	Ki67 Rabbit pAb
Remarks/Alias	KiA; MIB-; MIB-1; PPP1R105; Ki67
Species	Human
GeneID	4288
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 460-680 of human Ki67 (NP_002408.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P46013
Protein Weight	359 kDa
Shipping	Ice bag