

Anti-Kiaa0513 Rabbit Polyclonal Antibody For Wb, Elisa - O60268

Item Number: CM0019122

Introduction

Rabbit polyclonal antibody against human KIAA0513 (Uncharacterized protein KIAA0513), validated for Western blot and ELISA. Cross-reacts with mouse and rat, suitable for neuroscience and brain protein studies.

[Learn More](#)

Parameter	Specification
Product Name	KIAA0513 Rabbit pAb
Remarks/Alias	KIAA0513
Species	Human
Gene ID (Human)	9764
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-301 of human KIAA0513 (NP_001273494.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O60268
Protein Weight	47kDa
Shipping	Ice bag