

Anti-Kif22 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q14807

Item Number: CM0023540



Introduction

This rabbit polyclonal antibody against KIF22 (KID, KNSL4) detects human and mouse samples. Suitable for Western blot, immunohistochemistry, and ELISA. Targets a mitotic kinesin-like protein; useful for cell division and cancer research.

[Learn More](#)

Parameter	Value
Product Name	KIF22 Rabbit pAb
Remarks/Alias	KID; OBP; KNSL4; OBP-1; OBP-2; SEMDJL2; A-328A3.2; KIF22
Species	Human
GeneID (Human)	3835
GeneID	3835
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 456-665 of human KIF22 (NP_015556.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q14807
Protein Weight	73kDa
Shipping	Ice bag