


Anti-Kir2DL3 Polyclonal Antibody For Wb, If/Icc, Elisa - P43628

Item Number: CM0022141



Introduction

High-quality rabbit polyclonal antibody against human KIR2DL3, validated for WB, IF/ICC, and ELISA. Suitable for NK cell receptor studies, immune checkpoint research. Cross-reacts with human and mouse.

[Learn More](#)

Parameter	Value
Product Name	KIR2DL3 Rabbit pAb
Remarks/Alias	p58, NKAT, GL183, NKAT2, CD158b, KIR2DL, NKAT2A, NKAT2B, CD158B2, KIR-K7b, KIR-K7c, KIR2DS5, KIRCL23, KIR-023GB, KIR2DL3
Species	Human
GeneID (human)	3804
GeneID	3804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-245 of human KIR2DL3 (NP_056952.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P43628
Protein Weight	38kDa
Shipping	Ice bag